











SPECIMEN COLLECTION GUIDE - MOLECULAR

<p><u>THIN PREP PAP</u> ORDER#144598</p>	<p><u>Abbott Multi-Collect Specimen Collection Kit</u></p>
 <p>⇒ HPV ⇒ GC/Chlamydia</p>	 <p>ORDER# 176526 Intended for the collection and transportation of male and female, swab and urine specimens for the detection of Chlamydia trachomatis and Neisseria gonorrhoeae</p>
<p><u>E-SWAB</u> ⇒ MRSA by PCR</p>  <p>ORDER# 132433 ORDER# 193216</p>	<p><u>BD Affirm Vaginitis Panel</u> ORDER# 108401</p> <ul style="list-style-type: none"> ⇒ Gardnerella ⇒ Trichomonas ⇒ Candida 
<p><u>E-SWAB</u> ORDER# 132433 ⇒ Group B Strep by PCR Vaginal/Rectal eSwab</p> 	<p><u>STERILE SCREW CAP CUP (no additive)</u> ORDER# 153706</p> <ul style="list-style-type: none"> ⇒ Nasopharyngeal washes: Flu A / Flu B / RSV virus panel Flu A / H1N1 2009 virus panel B pertussis and parapertussis FilmArray 20 Pathogen Respiratory panel ⇒ Bronchial wash - HSV 1&2 ⇒ Mycoplasma pneumoniae Nasophayngeal wash ⇒ Stool: C. diff by PCR ⇒ CMV Qualitative PCR: CSF - Urine - Bronchial Wash 
<p>BD UTM (Universal Transport Media) ORDER# 111005 By PCR:</p> <ul style="list-style-type: none"> ⇒ Varicella virus ⇒ Herpes virus 1&2 ⇒ Flu A / Flu B / RSV virus panel ⇒ Flu A / H1N1 2009 virus panel ⇒ Mycoplasma pneumoniae ⇒ B pertussis and parapertussis ⇒ Respiratory panel 	<p><u>UroVysion</u> ORDER# 135884 ⇒ For bladder cancer test collection use Urocyte Kit</p> 
<p><u>PURPLE TOP TUBE - WHOLE BLOOD</u> ORDER# 9603</p> <ul style="list-style-type: none"> ⇒ Factor V Leiden Mutation by Real Time PCR ⇒ Factor II Mutation by Real Time PCR ⇒ Cystic Fibrosis screen: original tube only ⇒ CD4/CD8 (T-Cell Subsets) ⇒ Lymphoma/Leukemia Immunophenotyping ⇒ Fetomaternal Bleed ⇒ Flow Cytometry ⇒ HLA B27 Detection by Flowcytometry 	<p><u>GREEN TOP TUBE SODIUM HEPARIN</u> ORDER# 9604 ⇒ Hematological FISH testing for peripheral blood or bone marrow aspirate</p> 
<p><u>PURPLE TOP TUBE - PLASMA</u> ORDER# 9603</p> <ul style="list-style-type: none"> ⇒ CMV Qualitative PCR ⇒ CMV Quantitative PCR ⇒ HCV Quantitative by PCR ⇒ HIV Quantitative by PCR ⇒ HBV QT PCR Quantitative 