



Chlamydia/Gonorrhea Collection Device

Leaking Container Action Notification

PCL Alverno has transitioned to a new method for Molecular CTNG Testing. We are now using the Abbott® m2000™ system. One of the benefits of the Abbott system is the single collection kit for all sample types.

There have been instances of “leaking” unused multi-Collect sample collection kits. It is important to inspect the collection kits before using. This leaking may be obvious – or may simply appear as salt crystals around the rim of the collection tube. If either of these conditions exist, please do not use the collection kit. The supplier is actively working to correct the issue.

Please keep track of any leaking kits, as they will be replaced at no charge.

Following is an example of the appearance of a leaking tube:



Please do not hesitate to call PCL Alverno directly with any questions at 1.800.937.5521 with any questions.

Thank you for your support with this conversion.