







Test Summary
Utilization
GUIDELINES

GLUCOSE TOLERANCE TESTING

Effective January 15, 2019

PCL Alverno offers testing to screen for Type 2 diabetes in gestational and non-gestational patients. PCL Alverno will only offer testing according to guidelines recommended by the American Diabetes Association (ADA) and the American College of Obstetrics and Gynecology (ACOG).

For non-gestational screening, which includes non-pregnant women, men, and children, the ADA recommends a two-hour glucose tolerance test for selected patients with evidence of impaired glucose tolerance either on fasting or random measurements, or in the presence of elevated hemoglobin A1C.

For gestational screening of women between 24 – 28 weeks of pregnancy, the ACOG prefers the two-step approach, where an abnormal one-hour challenge with 50g of oral glucose solution is followed up with a three-hour oral glucose tolerance test (OGTT). Both the ADA and the ACOG list a one-step gestational screen that can be used as an alternative to screen pregnant women as well. The ADA has noted that measurement of hemoglobin A1C also can be used, but it may not be suitable for use alone because of decreased sensitivity compared with the OGTT approaches. The criteria and components for each test are listed below. All other glucose tolerance testing will be discontinued as of **January 15, 2019**.

	Test Name	Dose of Oral Glucose Solution	Test Components	Reference Range
Two-Step Gestational Diabetes Screen	1-hour Gestational Diabetes Screen	50g	1-hour glucose	70 – 139 mg/dL
	3-hour Gestational OGTT	100g	Fasting glucose 1-hour glucose 2-hour glucose 3-hour glucose	70 – 94 mg/dL 70 – 179 mg/dL 70 – 154 mg/dL 70 – 139 mg/dL
Non-Gestational Diabetes Screen	2-hour Non-gestational OGTT	75g	Fasting glucose 2-hour glucose	70 – 99 mg/dL 70 – 139 mg/dL
One-Step Gestational Diabetes Screen	2-hour Gestational OGTT	75g	Fasting glucose 1-hour glucose 2-hour glucose	70 – 91 mg/dL 70 – 179 mg/dL 70 – 152 mg/dL

FIND US ON FACEBOOK, LINKEDIN AND GOOGLE PLUS!

2434 Interstate Plaza Drive, Hammond, IN 46324 toll free 800.937.5521 · phone 219.989.3700 · fax 219.989.3900